


PIPETTE TIPS - "economy"

- Manufactured from high quality pure polypropylene which is highly resistant to chemicals and thermal effects.
- Produced to match wide variety of pipettors for secure sealing and easy tip ejecting.
- Uniformity in size and shape allows perfect fitting with almost all micropipettes. Flashfree and uniform delivery orifice for consistent accuracy and reproducible sample delivery.
- Thin wall and small surface area at the tip minimize the possible tip retaining samples droplets.

catalogue number	capacity	colour	pack quantity
005.01.001	10 μ l	clear	1000 pcs
005.01.002	200 μ l	yellow	1000 pcs
005.01.003	1000 μ l	blue	500 pcs
005.01.004	5 ml	clear	250 pcs
005.01.005	10 ml	clear	100 pcs



RESERVOIR - "2 in 1 - non sterile"

- Reservoir has a standard 50 ml basin on one side and 12 individual 5 ml basins on the other side.
- Wells are angled both downward and outward, guiding multichannel pipettor tips toward a centered position, which help to avoid contamination.
- Imprinted indices identify samples and loose-fitting lids prevent evaporation and contamination.

catalogue number	pack quantity
006.07.001	1 piece



RESERVOIR - "3 in 1 - non sterile"

- Manufactured from polypropylene as a single reservoir with 3 layouts as 2 sides with 3 compartments.
- Ideal for use when filling multi tube blocks, microwell plates etc since this two side reservoir is actually three reservoirs in one.
- One side has a single channel reservoir of 50 ml with clearly visible volume markings at 10, 25 and 40 ml levels.
- The other side has two separate compartments such as: 8 separate channels with 5 ml volume each division, 12 separate channels with 5 ml volume each division.

catalogue number	pack quantity
006.08.001	1 piece