

**SEALING TAPE - "for petri dishes"**

- An ideal sealing tape for petri dishes. Waterproof.
- Chemically inert, resists most organic solvents and caustic agents.
- Remains flexible in moderate freezer temperatures and withstands intermittent exposure to 100°C.
- Stretches, conforms and seals the irregular perimeter of petri dishes.

| catalogue number | tape width | pack quantity |
|------------------|------------|---------------|
| 058.07.001       | 13 mm      | 1 piece       |

**PETRI DISH STAND - "acrylic"**

- Petri dish stand is very useful during inoculation, incubation and storage.
- Visibility from all sides enable the culture to be observed during incubation.
- Safely transports, incubates and stores petri dishes. Also helps streamline benchtop sample processing.
- Center divider provides space for labelling and incorporates built in finger grip for easy transport.

| catalogue number | for petri dishes | dish capacity | quantity in pack |
|------------------|------------------|---------------|------------------|
| 081.03.060       | 60 mm diameter   | 56 pieces     | 1 piece          |
| 081.03.090       | 90 mm diameter   | 60 pieces     | 1 piece          |

**PETRI DISH BASKET - "st. steel"**

- Petri dish basket is manufactured from stainless steel.
- A very useful and practical holder for transportation, autoclaving process of petri dishes upto 100 mm diameter.

| catalogue number | dimensions         | pack quantity |
|------------------|--------------------|---------------|
| 081.20.001       | 114 x 102 x 245 mm | 1 piece       |

