



### PIPETTE PUMP - "electrical"

- Manufactured from UV resistant resin powered with a long lasting motor for quick but also sensitive aspiration and dispensing of liquids. Can be used with all pipettes from 0,1 ml to 100 ml volume.
- With the help of the powerfull, quiet and vibration-free motor a 25 ml pipette can be filled in less than 5 seconds.
- The ergonomic, light weight and well balanced handle offers comfort and fatigue-free pipetting.
- Hand-neutral design for right and left hand users, concave padded triggers, and convenient controls ensure the filler to perform many hours of stress-free filling and dispensing.
- 8 different level of aspiration and dispensing speed. Speed is easily set with a single switch.
- Large LCD display provides visual confirmation of the selected filling speed.
- Force liquid out of pipettes by setting the unit to low or high speed blow-out mode, or flip to gravity mode for use with "to contain" pipettes and dispense the liquid without motor.
- With the help of fully autoclavable silicone adapter pipettes are firmly hold in position.
- A replaceable 0,45 µl hydrophobic membrane filter in the intake cone secures contaminations and prevents the impurities to penetrate into the interior.
- Efficient lithium-ion battery offers long runtime such as 8-9 hours of continuous use.
- Batteries can be recharged within 2-3 hours while in use, or while resting in the desktop holder. With battery alarm.
- Large LCD display provides visual confirmation of remaining battery power.

catalogue number	for pipettes	aspirate speed	dispense speed	filter	batterie	batterie service life	charging time	pack quantity
010.01.005	0,1 - 100 ml	8 speed	8 speed	0,45 µm	lithium - ion	more than 8 hours	2 - 3 hours	1 piece