

**ADAPTERS - "with thread"**

- Manufactured from borosilicate glass 3.3.
- Excellent, smooth grinding of joints prevents airleaks.

**ADAPTERS - "with thread & cap"**

- Made from borosilicate glass 3.3 and they are ideal to be used with thermometers with joint glassware.
- Screw joint provides easy depth adjustment of thermometer.
- Can be used with thermometers with 6 - 8 mm diameter.



adapter straight	cone NS size	neck thread	pack quantity
065.19.141	NS 14/23	GL 14	1 piece
065.19.142	NS 14/23	GL 18	1 piece
065.19.291	NS 29/32	GL 14	1 piece
065.19.292	NS 29/32	GL 18	1 piece
065.19.293	NS 29/32	GL 25	1 piece

catalogue number	form	neck thread	cone NS size	pack quantity
065.21.014	straight	GL 14	NS 14/23	1 piece
065.21.029	straight	GL 14	NS 29/32	1 piece

**ADAPTERS - "receiver"**

- Made from borosilicate glass 3.3. Suitable for the usage of assembling various reactor systems.
- Available 2 types:
  - bend with angle 75°
  - bend with angle 90°

adapter bend-75°	adapter bend-90°	cone NS size	pack quantity
065.37.014	065.38.014	NS 14/23	1 piece
065.37.029	065.38.029	NS 29/32	1 piece

**ADAPTERS - "distillation"**

- Made from borosilicate glass 3.3. Suitable for feeding distillate from condensers into receiving flasks.
- Socket inclined 105° with an outer tube of 12 mm diameter.
- Available 2 types:
  - short stem - 65 mm.
  - long stem - 200 mm.

adapters 65 mm	adapters 200 mm	form	socket NS size	pack quantity
065.30.014	065.31.014	bend	NS 14/23	1 piece
065.30.029	065.31.029	bend	NS 29/32	1 piece

