



JOINT CLIPS - "POM"

- Manufactured from POM which has excellent chemical resistance as well as a very good mechanical strength.
- Perfect tools to keep two ground joints fastened firmly to each other.

catalogue number	POM colour	NS size	pack quantity
063.04.010	black	NS 10/19	1 piece
063.04.012	purple	NS 12/21	1 piece
063.04.014	yellow	NS 14/23	1 piece
063.04.019	blue	NS 19/26	1 piece
063.04.024	green	NS 24/29	1 piece
063.04.029	red	NS 29/32	1 piece
063.04.034	orange	NS 34/35	1 piece
063.04.045	brown	NS 45/40	1 piece



JOINT CLIPS - "for spherical joints"

- Manufactured from best stainless steel. Suitable only for KS 29 joints.
- A perfect tool to hold the spherical joints firmly during laboratory applications.
- Adjustable locking screw enables holding the clips at desired position.

catalogue number	material	pack quantity
063.05.001	stainless steel	1 piece



JOINT GREASE - "for stopcock & joints"

- High vacuum quality joint grease without silicone water washable and water soluble, non solidifying.
- Ideal for sealing and greasing or lubricating ground and spherical joints, plastic and metal connectors, stopcocks etc.
- Solves the most demanding lubricating problems in precision engineering and industry.
- Only a very small amount as very thin layer application is required to achieve the best result.
- Supplied in tubes of 60 gram lubricant weight.

catalogue number	product description	color	working temperature	pack quantity
063.16.001	normal vacuum	clear	-40° C to 300° C	1 piece
063.16.002	high vacuum	brown	-40° C to 320° C	1 piece