

BEAKERS - "PTFE" - "stainless steel" & "quartz glass"**"PTFE" - "low form"**


- Manufactured from high quality PTFE with excellent resistance to all chemicals.
- Manufactured from molded, pure PTFE with a smooth internal finish.
- Flat base enables a flat contact to the surface for effective heat transfer.
- PTFE beakers have excellent temperature stability from -200°C to +260°C.

catalogue number	volume	pack quantity
025.03.050	50 ml	1 piece
025.03.100	100 ml	1 piece
025.03.250	250 ml	1 piece
025.03.500	500 ml	1 piece
025.03.901	1000 ml	1 piece

**"st. steel" - "low form"**

- Manufactured from high quality AISI 304 stainless steel for improved corrosion resistance.
- Robust seamless construction for long working life.
- All surfaces are smooth and being highly polished enabling easy cleaning.
- Suitable for pharmaceutical laboratory applications since all surfaces are non-reacting with chemicals.

catalogue number	volume	pack quantity
025.04.100	100 ml	1 piece
025.04.250	250 ml	1 piece
025.04.500	500 ml	1 piece
025.04.901	1000 ml	1 piece

**"quartz glass" - "low form"**

- Manufactured from pure quartz and has a clear appearance like a glass.
- Designed to hold and pour chemical substances even in high temperatures up to 1050°C.
- Reinforced rims and spouts to prevent breakage during washing and transportation.

catalogue number	volume	pack quantity
025.06.010	10 ml	1 piece
025.06.025	25 ml	1 piece
025.06.050	50 ml	1 piece
025.06.100	100 ml	1 piece
025.06.250	250 ml	1 piece