

BEAKERS WITH HANDLE - "polypropylene"

*"embossed scale"**"blue printed scale"*

➤ Manufactured from ultra clear virgin polypropylene which provides thick walls with superior chemical resistance according to ISO 7056.

➤ They have specially designed functional spouts in order to minimize possible spills.

➤ Due to the unique and ergonomic shape of handles these beakers are stackable for space saving purpose.

➤ With embossed volume graduations and inscriptions.

➤ Manufactured from ultra clear virgin polypropylene which provides thick walls with superior chemical resistance according to ISO 7056.

➤ They have specially designed functional spouts in order to minimize possible spills.

➤ Due to the unique and ergonomic shape of handles these beakers are stackable for space saving purpose.

➤ With blue printed volume graduations and inscriptions.

embossed scale	volume	pack quantity
026.02.250	250 ml	1 piece
026.02.500	500 ml	1 piece
026.02.901	1000 ml	1 piece
026.02.902	2000 ml	1 piece
026.02.903	3000 ml	1 piece

blue printed scale	volume	pack quantity
026.04.250	250 ml	1 piece
026.04.500	500 ml	1 piece
026.04.901	1000 ml	1 piece
026.04.902	2000 ml	1 piece
026.04.903	3000 ml	1 piece



Due to the unique and ergonomic shape of handles these beakers are stackable for space saving purpose.



Blue printed volume graduations and inscriptions provides an excellent view.