


## ERLENMEYERS - "glass" - "narrow neck" &amp; "ground neck"



"clear" - "narrow neck"

- Manufactured from borosilicate glass 3.3. According to ISO 1773.
- Reinforced rims to prevent breaking during washing and transportation.
- Inscriptions are in white enamel.


catalogue number	volume	pack quantity
027.01.025	25 ml	10 pieces
027.01.050	50 ml	10 pieces
027.01.100	100 ml	10 pieces
027.01.250	250 ml	10 pieces
027.01.300	300 ml	10 pieces
027.01.500	500 ml	10 pieces
027.01.901	1000 ml	10 pieces
027.01.902	2000 ml	6 pieces
027.01.903	3000 ml	1 piece
027.01.905	5000 ml	1 piece



"amber" - "ground neck"

- Manufactured from borosilicate glass 3.3. According to ISO 1773.
- Reinforced rims to prevent breaking during washing and transportation.
- Inscriptions are in white enamel.

catalogue number	volume	neck NS size	pack quantity
028.11.053	50 ml	NS 29/32	1 piece
028.11.103	100 ml	NS 29/32	1 piece
028.11.252	250 ml	NS 29/32	1 piece
028.11.502	500 ml	NS 29/32	1 piece
028.11.901	1000 ml	NS 29/32	1 piece
028.11.902	2000 ml	NS 29/32	1 piece



"clear" - "ground neck"

- Manufactured from borosilicate glass 3.3. According to ISO 1773.
- They meet the highest quality standards according to DIN EN ISO 4797.
- NS necks are made in compliance with DIN 12242 standards.
- Inscriptions are in white enamel

catalogue number	volume	neck NS size	pack quantity
028.01.050	50 ml	NS 14/23	1 piece
028.01.051	50 ml	NS 19/26	1 piece
028.01.052	50 ml	NS 24/29	1 piece
028.01.053	50 ml	NS 29/32	1 piece
028.01.100	100 ml	NS 14/23	1 piece
028.01.101	100 ml	NS 19/26	1 piece
028.01.102	100 ml	NS 24/29	1 piece
028.01.103	100 ml	NS 29/32	1 piece
028.01.250	250 ml	NS 19/26	1 piece
028.01.251	250 ml	NS 24/29	1 piece
028.01.252	250 ml	NS 29/32	1 piece
028.01.300	300 ml	NS 29/32	1 piece
028.01.500	500 ml	NS 19/26	1 piece
028.01.501	500 ml	NS 24/29	1 piece
028.01.502	500 ml	NS 29/32	1 piece
028.01.901	1000 ml	NS 29/32	1 piece
028.01.902	2000 ml	NS 29/32	1 piece



Glass stoppers for ground neck erlenmeyers can be visited at page 110.

Aluminium covers for erlenmeyers can be visited at page 110.