

VOLUMETRIC FLASKS - "polypropylene" - "class B"

- Manufactured from glass-like transparent special grade polypropylene.
- They offer thick walls which eliminates distortion, high transparency, chemical resistance and remarkable heat tolerance.
- Calibrated to contain (In) and the tolerances are according to DIN/ISO Class "B" standards.
- They meet the highest quality standards according to ISO 1042.
- Flasks have inscriptions and meniscus ring mark in blue enamel which is specially designed for plastic and resistant to 60°C.
- For consistent accuracy, P.P volumetric flasks should not be used with liquids over 60°C and or kept under direct sunlight or overheated conditions.



clear flasks blue scale	volume	tolerance	approximate height	neck inner diameter	neck size	pack quantity
014.03.025	25 ml	± 0.080 ml	110 mm	10.00 mm	NS 10/19	1 piece
014.03.050	50 ml	± 0.120 ml	140 mm	12.50 mm	NS 12/21	1 piece
014.03.100	100 ml	± 0.200 ml	170 mm	14.50 mm	NS 14/23	1 piece
014.03.250	250 ml	± 0.300 ml	220 mm	14.50 mm	NS 14/23	1 piece
014.03.500	500 ml	± 0.500 ml	260 mm	18.80 mm	NS 19/26	1 piece
014.03.901	1000 ml	± 0.800 ml	300 mm	24.00 mm	NS 24/29	1 piece