



FLASKS - "flat & round bottom" - "glass"

- Manufactured from borosilicate glass 3.3. According to ISO 1773 and DIN 12353.
- Reinforced rims are used to prevent breakage during washing and transportation.

flat bottom flasks	round bottom flasks	volume	approx. diameter	pack quantity
029.01.050	029.02.050	50 ml	51 mm	10 pieces
029.01.100	029.02.100	100 ml	64 mm	10 pieces
029.01.250	029.02.250	250 ml	85 mm	10 pieces
029.01.500	029.02.500	500 ml	105 mm	10 pieces
029.01.901	029.02.901	1000 ml	131 mm	10 pieces
029.01.902	029.02.902	2000 ml	166 mm	6 pieces
029.01.904	029.02.904	4000 ml	207 mm	1 piece
029.01.906	029.02.906	6000 ml	236 mm	1 piece
029.01.910	029.02.910	10000 ml	279 mm	1 piece

FLASKS - "NS neck"

- Made from borosilicate glass 3.3. They meet the highest quality standards according to DIN 12353.
- NS necks are made in compliance with DIN 12242 standards. Inscriptions are in white enamel which is specially designed for glass.

flasks (NS) flat bottom	flasks (NS) round bottom	volume	neck NS size	pack quantity
030.01.050	030.02.050	50 ml	NS 14/23	1 piece
030.01.051	030.02.051	50 ml	NS 19/26	1 piece
030.01.052	030.02.052	50 ml	NS 24/29	1 piece
030.01.053	030.02.053	50 ml	NS 29/32	1 piece
030.01.100	030.02.100	100 ml	NS 14/23	1 piece
030.01.101	030.02.101	100 ml	NS 19/26	1 piece
030.01.102	030.02.102	100 ml	NS 24/29	1 piece
030.01.103	030.02.103	100 ml	NS 29/32	1 piece
030.01.250	030.02.250	250 ml	NS 19/26	1 piece
030.01.251	030.02.251	250 ml	NS 24/29	1 piece
030.01.252	030.02.252	250 ml	NS 29/32	1 piece
030.01.500	030.02.500	500 ml	NS 24/29	1 piece
030.01.501	030.02.501	500 ml	NS 29/32	1 piece
030.01.900	030.02.900	1000 ml	NS 24/29	1 piece
030.01.901	030.02.901	1000 ml	NS 29/32	1 piece
030.01.902	030.02.902	2000 ml	NS 29/32	1 piece
030.01.904	030.02.904	4000 ml	NS 29/32	1 piece



FLASKS - "round bottom" - "NS neck" - "amber"

- Made from amber colour borosilicate glass 3.3. They meet the highest quality standards according to DIN 12348.
- NS necks are made in compliance with DIN 12242 standards. Inscriptions are in white enamel which is specially designed for glass.

catalogue number	volume	neck NS size	pack quantity
030.12.053	50 ml	NS 29/32	1 piece
030.12.103	100 ml	NS 29/32	1 piece
030.12.252	250 ml	NS 29/32	1 piece
030.12.501	500 ml	NS 29/32	1 piece
030.12.901	1000 ml	NS 29/32	1 piece