



FLASKS - "ground neck" - "2 neck"

➤ Manufactured from borosilicate glass 3.3 according to DIN 12394 with thick walls, heavy bases and reinforced rims which prevent breakage.

➤ Available 2 types: ➤ parallel side neck
➤ angled side neck

angled side arm	parallel side arm	volume	neck NS size	side NS size
030.24.050	030.30.050	50 ml	NS 14/23	NS 14/23
030.24.100	030.30.100	100 ml	NS 14/23	NS 14/23
030.24.250	030.30.250	250 ml	NS 14/23	NS 14/23
030.25.100	030.31.100	100 ml	NS 24/29	NS 14/23
030.25.250	030.31.250	250 ml	NS 24/29	NS 14/23
030.25.500	030.31.500	500 ml	NS 24/29	NS 14/23
030.25.901	030.31.901	1000 ml	NS 24/29	NS 14/23
030.26.100	030.32.100	100 ml	NS 24/29	NS 19/26
030.26.250	030.32.250	250 ml	NS 24/29	NS 19/26
030.26.500	030.32.500	500 ml	NS 24/29	NS 19/26
030.26.901	030.32.901	1000 ml	NS 24/29	NS 19/26
030.27.100	030.33.100	100 ml	NS 29/32	NS 14/23
030.27.250	030.33.250	250 ml	NS 29/32	NS 14/23
030.27.500	030.33.500	500 ml	NS 29/32	NS 14/23
030.27.901	030.33.901	1000 ml	NS 29/32	NS 14/23
030.27.902	030.33.902	2000 ml	NS 29/32	NS 14/23
030.28.100	030.34.100	100 ml	NS 29/32	NS 19/26
030.28.250	030.34.250	250 ml	NS 29/32	NS 19/26
030.28.500	030.34.500	500 ml	NS 29/32	NS 19/26
030.28.901	030.34.901	1000 ml	NS 29/32	NS 19/26
030.28.902	030.34.902	2000 ml	NS 29/32	NS 19/26
030.29.100	030.35.100	100 ml	NS 29/32	NS 29/32
030.29.250	030.35.250	250 ml	NS 29/32	NS 29/32
030.29.500	030.35.500	500 ml	NS 29/32	NS 29/32
030.29.901	030.35.901	1000 ml	NS 29/32	NS 29/32
030.29.902	030.35.902	2000 ml	NS 29/32	NS 29/32



FLASKS - "ground neck" - "3 neck"

➤ Manufactured from borosilicate glass 3.3 according to DIN 12394 with thick walls, heavy bases and reinforced rims which prevent breakage.

➤ Available 2 types: ➤ parallel side neck
➤ angled side neck

angled side arm	parallel side arm	volume	neck NS size	side NS size
030.36.250	030.42.250	250 ml	NS 14/23	NS 14/23
030.37.250	030.43.250	250 ml	NS 24/29	NS 14/23
030.37.500	030.43.500	500 ml	NS 24/29	NS 14/23
030.37.901	030.43.901	1000 ml	NS 24/29	NS 14/23
030.38.250	030.44.250	250 ml	NS 24/29	NS 19/26
030.38.500	030.44.500	500 ml	NS 24/29	NS 19/26
030.38.901	030.44.901	1000 ml	NS 24/29	NS 19/26
030.38.902	030.44.902	2000 ml	NS 24/29	NS 19/26
030.39.250	030.45.250	250 ml	NS 29/32	NS 14/23
030.39.500	030.45.500	500 ml	NS 29/32	NS 14/23
030.39.901	030.45.901	1000 ml	NS 29/32	NS 14/23
030.40.250	030.46.250	250 ml	NS 29/32	NS 19/26
030.40.500	030.46.500	500 ml	NS 29/32	NS 19/26
030.40.901	030.46.901	1000 ml	NS 29/32	NS 19/26
030.41.250	030.47.250	250 ml	NS 29/32	NS 29/32
030.41.500	030.47.500	500 ml	NS 29/32	NS 29/32
030.41.901	030.47.901	1000 ml	NS 29/32	NS 29/32
030.41.902	030.47.902	2000 ml	NS 29/32	NS 29/32