



### FLASKS - "kjeldahl"

➤ Made from borosilicate glass 3.3 and meeting quality standards according to DIN 12360.

catalogue number	neck NS size	volume	pack quantity
030.57.100	NS 19/26	100 ml	1 piece
030.57.250	NS 29/32	250 ml	1 piece
030.57.500	NS 29/32	500 ml	1 piece
030.57.750	NS 29/32	750 ml	1 piece
030.57.901	NS 29/32	1000 ml	1 piece

### FLASKS - "nitrogen"

➤ Manufactured from borosilicate glass 3.3 and meeting quality standards according to DIN and ISO standards.

➤ NS necks are in compliance with DIN 12242 standards. They have angled side arm with PTFE stopcock.

catalogue number	nitrogen volume	neck NS size	pack quantity
030.58.025	25 ml	NS 14/23	1 piece
030.58.050	50 ml	NS 14/23	1 piece
030.58.100	100 ml	NS 14/23	1 piece
030.58.250	250 ml	NS 14/23	1 piece
030.59.050	50 ml	NS 29/32	1 piece
030.59.100	100 ml	NS 29/32	1 piece
030.59.250	250 ml	NS 29/32	1 piece
030.59.500	500 ml	NS 29/32	1 piece
030.59.901	1000 ml	NS 29/32	1 piece



### FLASK TONG



➤ Manufactured from chromated steel. Insulated jaws are covered with special coating with perfect insulation enabling to hold all types of vessels in most secure way.

catalogue number	total length	holding capacity	pack quantity
029.03.001	290 ± 10 mm	1000 ml	1 piece

### FLASK TONG



➤ Manufactured from chromated steel. Insulated jaws are covered with special coating with perfect insulation enabling to hold all types of vessels in most secure way.

catalogue number	total length	holding capacity	pack quantity
029.03.002	280 ± 10 mm	1000 ml	1 piece



### FLASK STANDS - "cork" & "polypropylene"

➤ Cork flask stands are manufactured from the best quality high density cork.

➤ They are ideal tools to be used with all volumes of round bottom flasks.

catalogue number	recommended flask volume
029.04.001	50 - 100 ml
029.04.002	250 - 500 ml
029.04.003	1000 - 2000 ml
029.04.004	3000 ml
029.04.005	4000 ml

➤ Plastic flask stands are made from P.P.

➤ Unique wall design with stepped concentric rings offers an excellent support at every point of contact.

➤ Stands can be stacked together for convenient storage. Supplied 1 piece.

catalogue number	for flasks volume
029.14.902	100 - 2000 ml

