




CONDENSERS - "liebig"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

| with glass side arm | with P.P side arm | effective length | joint NS size | pack quantity |
|---------------------|-------------------|------------------|---------------|---------------|
| 033.11.160 | 033.12.160 | 160 mm | NS 14/23 | 1 piece |
| 033.11.250 | 033.12.250 | 250 mm | NS 14/23 | 1 piece |
| 033.11.252 | 033.12.252 | 250 mm | NS 24/29 | 1 piece |
| 033.11.253 | 033.12.253 | 250 mm | NS 29/32 | 1 piece |
| 033.11.400 | 033.12.400 | 400 mm | NS 14/23 | 1 piece |
| 033.11.402 | 033.12.402 | 400 mm | NS 24/29 | 1 piece |
| 033.11.403 | 033.12.403 | 400 mm | NS 29/32 | 1 piece |



CONDENSERS - "allihn"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

| with glass side arm | with P.P side arm | effective length | joint NS size | pack quantity |
|---------------------|-------------------|------------------|---------------|---------------|
| 033.13.160 | 033.14.160 | 160 mm | NS 14/23 | 1 piece |
| 033.13.250 | 033.14.250 | 250 mm | NS 14/23 | 1 piece |
| 033.13.252 | 033.14.252 | 250 mm | NS 24/29 | 1 piece |
| 033.13.253 | 033.14.253 | 250 mm | NS 29/32 | 1 piece |
| 033.13.400 | 033.14.400 | 400 mm | NS 14/23 | 1 piece |
| 033.13.402 | 033.14.402 | 400 mm | NS 24/29 | 1 piece |
| 033.13.403 | 033.14.403 | 400 mm | NS 29/32 | 1 piece |

