

## CONDENSERS - "coiled"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

with glass side arm	with P.P side arm	effective length	joint NS size	pack quantity
033.15.160	033.16.160	160 mm	NS 14/23	1 piece
033.15.253	033.16.253	250 mm	NS 29/32	1 piece
033.15.403	033.16.403	400 mm	NS 29/32	1 piece



## CONDENSERS - "dimroth"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

with glass side arm	with P.P side arm	effective length	joint NS size	pack quantity
033.17.160	033.18.160	160 mm	NS 14/23	1 piece
033.17.163	033.18.163	160 mm	NS 29/32	1 piece
033.17.250	033.18.250	250 mm	NS 14/23	1 piece
033.17.253	033.18.253	250 mm	NS 29/32	1 piece
033.17.400	033.18.400	400 mm	NS 14/23	1 piece
033.17.403	033.18.403	400 mm	NS 29/32	1 piece

