


## COLUMNS - "chromatography"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- NS necks are in compliance with DIN 12242 standards.
- Available 2 types:
  - without frit
  - with frit

column without frit	column with frit	effective length	inside diameter	socket NS size	pack quantity
066.10.001	066.11.001	200 mm	10 mm	NS 14/23	1 piece
	066.11.003	200 mm	15 mm	NS 14/23	1 piece
	066.11.002	300 mm	10 mm	NS 14/23	1 piece
066.10.003	066.11.004	400 mm	20 mm	NS 29/32	1 piece
066.10.004	066.11.005	600 mm	30 mm	NS 29/32	1 piece



Chromatographic columns with sintered frit, pore size 1.



## SOXHLET CONDENSERS - "allihn"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

with glass side arm	with P.P side arm	for extractor volume	cone NS size	pack quantity
035.11.030	035.12.030	30 ml	NS 29/32	1 piece
035.11.070	035.12.070	70 ml	NS 34/35	1 piece
035.11.250	035.12.250	100 / 250 ml	NS 45/40	1 piece
035.11.500	035.12.500	500 ml	NS 45/40	1 piece

