

## SOXHLET CONDENSERS - "dimroth"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They are carefully produced in order to meet the highest quality standards.
- Supplied with either glass or polypropylene side arm.
- Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing.
- The side arms accept tubings with 8/9 mm bore size.

with glass side arm	with P.P side arm	for extractor volume	cone NS size	pack quantity
035.13.030	035.14.030	30 ml	NS 29/32	1 piece
035.13.070	035.14.070	70 ml	NS 34/35	1 piece
035.13.250	035.14.250	100 / 250 ml	NS 45/40	1 piece
035.13.500	035.14.500	500 ml	NS 45/40	1 piece



## SOXHLET EXTRACTORS - "W/O stopcock"

- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- They meet the highest quality standards according to DIN 12602.

catalogue number	volume	pack quantity
036.01.030	30 ml	1 piece
036.01.070	70 ml	1 piece
036.01.100	100 ml	1 piece
036.01.150	150 ml	1 piece
036.01.250	250 ml	1 piece
036.01.500	500 ml	1 piece

catalogue number	socket NS size
036.01.030	NS 29/32
036.01.070	NS 34/35
036.01.100	NS 45/40
036.01.150	NS 45/40
036.01.250	NS 45/40
036.01.500	NS 45/40



## CELLULOSE THIMBLES

- Manufactured from carefully selected pure cellulose.
- Ideal for soxhlet or other extractions to obtain reproducible results without any harm to the samples.
- They have perfectly controlled dimensional accuracy which provides excellent fitting to the extractor.
- They are also suitable for extracting solid or liquid particles from the air or gasses.



catalogue number	diameter	height	pack quantity
102.28.080	28 mm	80 mm	25 pieces
102.30.080	30 mm	80 mm	25 pieces
102.33.080	33 mm	80 mm	25 pieces
102.33.094	33 mm	94 mm	25 pieces
102.34.120	34 mm	120 mm	25 pieces