



FUNNELS - "stainless steel"

- Manufactured from best quality polished stainless steel.
- Stem is vented to allow undisturbed flow of liquids.

catalogue number	funnel diameter	height	pack quantity
041.06.080	80 ± 5 mm	115 ± 5 mm	1 piece
041.06.100	100 ± 5 mm	125 ± 5 mm	1 piece
041.06.120	120 ± 5 mm	145 ± 5 mm	1 piece
041.06.150	150 ± 5 mm	190 ± 5 mm	1 piece



FUNNELS - "polyethylene" - "industrial"

- Manufactured from clear, virgin HD-PE.
- They are used to fill in drums, carboys and so on.

catalogue number	funnel diameter	pack quantity
041.02.200	200 mm	1 piece
041.02.250	250 mm	1 piece
041.02.300	300 mm	1 piece



FUNNELS - "for powder" - "polypropylene"

- Manufactured from clear, polypropylene and specially designed for powders. Can also be used with liquids.
- Large bore and parallel stem avoid bridging of powders during transfer.

catalogue number	funnel diameter	stem diameter	pack quantity
041.07.065	65 mm	15 mm	1 piece
041.07.080	80 mm	21 mm	1 piece
041.07.100	100 mm	22 mm	1 piece
041.07.120	120 mm	26 mm	1 piece



FUNNELS - "for powder" - "glass"

- Manufactured from borosilicate glass 3.3 with short stem and 60° cone angle.
- These funnels are practical and suitable for filling liquids and powders.

catalogue number	funnel diameter	cone NS size	pack quantity
041.08.060	60 mm	NS 14/23	1 piece
041.08.080	80 mm	NS 19/26	1 piece
041.08.100	100 mm	NS 24/29	1 piece
041.08.120	120 mm	NS 29/32	1 piece