



### FUNNEL HOLDERS - "stainless steel"

- Manufactured from best quality stainless steel.
- Unique ring shape with slit offers practical use with separating funnels.
- Suitable to hold all types of funnels with central diameter up to 150 mm.
- Available 2 types:
  - without boss head
  - with boss head

funnel holder W/O boss head	funnel holder with boss head	holder diameter	pack quantity
042.03.040	042.04.040	40 mm	1 piece
042.03.060	042.04.060	60 mm	1 piece
042.03.080	042.04.080	80 mm	1 piece
042.03.100	042.04.100	100 mm	1 piece



### FUNNEL STAND

- Made from stainless steel, and supplied with screw head to adjust the height of the holder.
- Adjustable holding system enables flexible and practical use of the funnel stand.
- Thanks to the adjustable height of holding part which offers practical use almost all sizes of funnels.

catalogue number	pack quantity
042.06.010	1 piece



Funnel holder without slit consists two pieces.



Funnel holder with slit consists single piece.

### FUNNEL HOLDERS - "polypropylene"

- Manufactured from polypropylene. Suitable to hold all types of laboratory funnels with diameter from 60 mm up to 150 mm.
- Two tapered wedges are used to ensure firm grip on rods having diameter 10 -13 mm.
- Available 2 types:
  - without slit
  - with slit
- Funnel holder with slit is mostly recommended for use with separating funnel.

funnel holder without slit	funnel holder with slit	pack quantity
042.05.001	042.06.001	1 piece