

FILTER FLASKS - "glass" & "polyethylene"



"glass" - "glass side arm"

- Manufactured from borosilicate glass 3.3 with thick walls according to ISO 6556 and DIN 12476 standards.
- Side arms are suitable to fit 6-10 mm hose diameter.

catalogue number	volume	pack quantity
040.01.100	100 ml	1 piece
040.01.250	250 ml	1 piece
040.01.500	500 ml	1 piece
040.01.901	1000 ml	1 piece
040.01.902	2000 ml	1 piece



"glass" - "P.P side arm"

- Manufactured from borosilicate glass 3.3 with thick walls according to ISO 6556 and DIN 12476 standards.
- Polypropylene side arm has a screw thread connector for easy and safe fixing of the tubings. Side arms are suitable to fit 6-10 mm hose diameter.

catalogue number	volume	pack quantity
040.02.100	100 ml	1 piece
040.02.250	250 ml	1 piece
040.02.500	500 ml	1 piece
040.02.901	1000 ml	1 piece
040.02.902	2000 ml	1 piece



"polyethylene"

- Manufactured from durable thick walled polyethylene.
- Side arms are suitable to fit 6-10 mm hose diameter.

catalogue number	volume	pack quantity
040.03.901	1000 ml	1 piece



GASKETS - "Guko"

- Manufactured from naturel rubber.
- Conical body form enables perfect fitting of filter funnels and crucibles in to suction flasks.
- The set consists of 7 pieces of gaskets suitable for all sizes of filter funnels.

catalogue number	pack quantity
040.06.001	1 set