


## VIALS & CAPS & SEPTAS



➤ Vials are a critical link to optimal analysis. Each septum is specified to complement the overall system and enhance chemical performance.

➤ ISOLAB vial and caps are specially formulated and constructed to provide optimum system performance with minimal coring and superior chemical inertness.

➤ The softness and hardness (durometer) is matched to the instrument and designed to work excellently well with syringes and samples.

➤ The size of the septum is critical to proper sealing. All ISOLAB septas are specifically sized for proper sealing.

➤ Center hole pad is made from high grade virgin silicone. Vial screw caps are offered either with a slit or without slit. ISOLAB caps perfectly fits to the neck of the vial.

➤ Manufactured from high quality polypropylene to exact manufacturing tolerances and lined in a controlled manufacturing environment, tested by head space GC to ensure proper curing of silicone polymer.

### SCREW NECK VIALS & CAPS - "N8" & "N9"



### CRIMP NECK VIALS & CAPS - "N11"



➤ Screw neck vials are ideal to be used for modern auto samplers.

➤ Volume graduation enables filling the sample into the vial to a desired level. Large marking area enables better sample identification.

➤ Crimp neck vials are ideal to be used for modern auto samplers.

➤ Volume graduation enables filling the sample into the vial to a desired level. Large marking area enables better sample identification.

screw neck vials	volume	neck N size	colour	vial dimensions	pack quantity
095.00.001	2.0 ml	N8	clear	11.6 × 32 mm	100 pieces
095.00.002	2.0 ml	N8	amber	11.6 × 32 mm	100 pieces
095.01.001	2.0 ml	N9	clear	11.6 × 32 mm	100 pieces
095.01.002	2.0 ml	N9	amber	11.6 × 32 mm	100 pieces

crimp neck vials	volume	neck N size	colour	vial dimensions	pack quantity
095.02.001	2.0 ml	N11	clear	11.6 × 32 mm	100 pieces
095.02.002	2.0 ml	N11	amber	11.6 × 32 mm	100 pieces

screw caps for vials	neck N size	septum description	pack quantity
096.00.001	N8	red silicone (w/o slit) / white PTFE	100 pieces
096.00.002	N8	red silicone (with slit) / white PTFE	100 pieces
096.01.001	N9	white silicone (w/o slit) / red PTFE	100 pieces
096.01.002	N9	white silicone (with slit) / blue PTFE	100 pieces

crimp caps for vials	neck N size	septum description	pack quantity
096.02.001	N11	red silicone (w/o slit) / white PTFE	100 pieces