

## STORAGE VIALS &amp; CAPS



storage vials	ml	neck N size	colour	vial dimensions	pack quantity
095.03.001	22 ml	N20	clear	23 x 85 mm	100 pieces
095.03.002	22 ml	N20	amber	23 x 85 mm	100 pieces
095.03.011	40 ml	N24	clear	28 x 95 mm	100 pieces
095.03.012	40 ml	N24	amber	28 x 95 mm	100 pieces

screw caps for vials	neck N size	septum description	pack quantity
096.05.001	N20	clear silicone (w/o hole) / white PTFE	100 pieces
096.05.011	N24	clear silicone (w/o hole) / white PTFE	100 pieces

## HEAD SPACE VIALS &amp; CAPS



head space vials	ml	neck N size	colour	vial dimensions	pack quantity
097.03.001	10 ml	N20	clear	22 x 46 mm	100 pieces
097.03.002	10 ml	N20	amber	22 x 46 mm	100 pieces
097.04.001	20 ml	N20	clear	22 x 75 mm	100 pieces
097.04.002	20 ml	N20	amber	22 x 75 mm	100 pieces

crimp caps for vials	neck N size	septum description	pack quantity
096.03.001	N20	silicone (w/o slit) / PTFE	100 pieces

## VIAL INSERTS



► Manufactured from clear glass with flat bottom having a volume of 0,25 ml or 0,30 ml.

► Vials inserts are a practical solution in applications where there is a limited sample volume.

► Also allows maximum sample recovery when used in conjunction with autosampler vials.

catalogue number	for vial size	dimensions (dia. x height)	volume	colour	pack quantity
097.05.108	N 8	5,0 x 31,5 mm	0,25 ml	clear	100 pieces
097.05.109	N 9	5,8 x 31,5 mm	0,30 ml	clear	100 pieces

## CRIMPER



catalogue no: 097.01.011

## DECAPPER



catalogue no: 097.02.011