

## ISO BOTTLES - "polypropylene"



➤ These impact resistant, easily recyclable, leak proof and excellently thermal stable bottles are autoclavable at max. 121°C.

➤ Manufactured from glass-like transparent special grade polypropylene.

➤ Angled neck form allows the liquid to flow out without any restrictions.

➤ GL 45 bottle neck provides easy filling and pouring. Unique neck design with moulded pouring ring prevents drops and liquid splash.

➤ Cap is made from polypropylene large ridges for increased grip with self sealing lock and integral security ring.

➤ Unique cap inner design ensures hermetical closure and guarantees an excellent leak and contamination proof bottle.

➤ Embossed volume graduation on side wall enables viewing of liquid level and reproducible filling of the bottle without use of cylinders or beakers.

➤ Sterile bottles are produced under germ-free sterile conditions with extreme protocols or sterility control then sent to gamma irradiation for further sterilisation.

clear bottle non-sterile	clear bottle sterile	bottle volume	pack quantity
061.15.100	061.16.100	100 ml	100 pieces / pack
061.15.250	061.16.250	250 ml	100 pieces / pack
061.15.500	061.16.500	500 ml	100 pieces / pack
061.15.901	061.16.901	1000 ml	50 pieces / pack
061.15.902	061.16.902	2000 ml	50 pieces / pack



Caps are made from robust and thick polypropylene with large ridges for increased grip.



Unique neck design with moulded pouring ring prevents drops and liquid splash.



Unique design of the cap with intended top plate guarantees tight and leakproof closure.



Integral security ring of the cap ensures safe and untouched protection of the bottle content.