

BOTTLES - "neutral glass"



- Manufactured from neutral glass which is relatively resistant to heat and almost all chemicals.
- Autoclavable at 121°C and 1.1 bar for 30 minutes. Caps should not be closed during autoclaving process.
- The angled finishing allows the liquid to flow without any restrictions.
- Cap is made from polypropylene with large ridges for increased grip.
- GL 45 bottle neck provides easy filling and pouring. Pouring ring prevents drops and liquid splash
- Inscriptions are in white enamel which is specially designed for glass.

catalogue number	volume	pack quantity
061.31.100	100 ml	10 pieces
061.31.250	250 ml	10 pieces
061.31.500	500 ml	10 pieces
061.31.901	1000 ml	10 pieces

ISO BOTTLES - "borosilicate glass"



- Manufactured from borosilicate glass 3.3 which is resistant to heat and almost all chemicals.
- Autoclavable at 121°C and 1.1 bar for 30 minutes. Caps should not be closed during autoclaving process.
- GL 45 bottle neck provides easy filling and pouring. Pouring ring prevents drops and liquid splash.
- Cap is made from polypropylene with large ridges for increased grip.
- The angled finishing allows the liquid to flow without any restrictions.
- Amber bottles prevent from UV lights to penetrate inside the bottle.
- Inscriptions are in white enamel which is specially designed for glass.

volume	clear bottles	pack quantity	amber bottles	pack quantity
100 ml	061.01.100	10 pieces	061.02.100	1 piece
250 ml	061.01.250	10 pieces	061.02.250	1 piece
500 ml	061.01.500	10 pieces	061.02.500	1 piece
1000 ml	061.01.901	10 pieces	061.02.901	1 piece
2000 ml	061.01.902	4 pieces	061.02.902	1 piece
3000 ml	061.01.903	1 piece	061.02.903	1 piece
5000 ml	061.01.905	1 piece	061.02.905	1 piece
10000 ml	061.01.910	1 piece	061.02.910	1 piece
20000 ml	061.01.920	1 piece	061.02.920	1 piece



Cap is made from polypropylene large ridges for increased grip with self sealing lock.



GL 45 bottle neck provides easy filling and pouring. Pouring ring prevents drops and liquid splash

