

BOTTLES - "polyethylene" - "wide neck"

- Manufactured from polyethylene as impact resistant, easily recyclable, leak proof and excellently thermal stability.
- Wide neck of the bottle enables easy filling and pouring.
- Clear bottles are manufactured from glass-like transparent special grade polyethylene.
- Guarantees perfect closure by an insert plug.
- The space between cap and insert plug can be used for securing private and confidential notes.
- Both at the bottle neck and at the cap there are holes which enable insertion of security seal wire.
- Embossed volume graduation on side wall enables viewing of liquid level and reproducible filling of the bottle without use of cylinders and beakers.

catalogue number	bottle volume	dimensions (w x d x h)	pack quantity
061.07.050	50 ml	37 x 39 x 74 mm	50 pieces
061.07.100	100 ml	45 x 50 x 90 mm	50 pieces
061.07.250	250 ml	62 x 57 x 113 mm	50 pieces
061.07.500	500 ml	83 x 71 x 140 mm	50 pieces
061.07.901	1000 ml	104 x 82 x 171 mm	50 pieces
061.07.902	2000 ml	134 x 107 x 210 mm	50 pieces



SECURITY SEAL & SEAL CUTTER

- Security seal is used ensure the safety locking and durability of the product inside the bottles.
- 4 different colors are used to identify among the seals with identification numbers printed on the seal plate.



Security seal cutter is a perfect tool to cut security seal wires easily.

catalogue number	pack quantity
061.20.001	1 piece



The space between cap and insert plug can be used for private and confidential notes.



Insert the seal wire through the holes of both the cap and bottle itself. Then insert the wire into the hole of the seal as shown.



Last step is to wind the seal key clockwise until the seal wire gets tight enough to provide optimum security lock.