

BOTTLES - "dropping" - "polyethylene"



➤ Manufactured from low density polyethylene.

➤ Dispensing tip provides regular drops at all times as well as provides continuous flow with the fine jet.

➤ The tip of the delivery tube can be cut in order to increase the drop volume.

➤ Amber bottles prevents UV light to penetrate inside the bottle.

clear bottle	amber bottle	bottle volume	pack quantity
062.04.100	062.04.101	100 ml	1 piece



Dispensing tip provides regular drops at all times as well as provides continuous flow with the fine jet.



➤ Manufactured from a mixture of high and low density polyethylene.

➤ Dispensing tip provides regular drops at all times as well as provides continuous flow with the fine jet.

➤ The tip of the delivery tube can be cut in order to increase the drop volume.

➤ Supplied with a push open seal cap that enables reclosing after using the dropping bottle.

catalogue number	volume	dimensions Ø x height	pack quantity
062.10.010	10 ml	21 x 69 mm	10 pieces
062.10.015	15 ml	23 x 76 mm	10 pieces
062.10.020	20 ml	25 x 82 mm	10 pieces
062.10.030	30 ml	30 x 82 mm	10 pieces



Dispensing tip provides regular drops at all times as well as provides continuous flow with the fine jet.

BOTTLES - "spray"



➤ Manufactured from PET with a long-life hand pump.

➤ Liquids can be finely vaporised with the pull of the hand pump trigger.

➤ Media contacting parts are polypropylene, polyethylene and stainless steel.

➤ Safety lock at the hand pump trigger prevents accidental use of bottle.

catalogue number	volume	pack quantity
062.11.250	250 ml	1 piece
062.11.500	500 ml	1 piece



Extremely homogeneous vaporisation is guaranteed with the pull of the trigger.