

WASH BOTTLES - "narrow neck" - "wide neck" - "integral"

**"narrow neck"**

- Manufactured from a mixture of high and low density polyethylene as semirigid, translucent and with high clarity.
- Very easy to use by pressing gently. The "swan neck" delivery tip produces a fine, straight flow without lateral spatter.
- The flow rate can be easily controlled by hand pressure.

Bottles can be fully emptied due to the movable outlet tube.

- Dispensing tip provides extra fine stream. The tip of the delivery tube can be cut in order to increase the liquid flow rate.

catalogue number	volume	pack quantity
062.05.100	100 ml	1 piece
062.05.250	250 ml	1 piece
062.05.500	500 ml	1 piece
062.05.901	1000 ml	1 piece

**"wide neck"**

- Manufactured from a mixture of high and low density polyethylene as semi rigid, and with very high clarity.
- Wide neck for easy filling of the liquid. Colour coded caps enable easy identification of the bottle content.
- Ribbed cap design allows easy opening and closing. One piece cap with angled delivery tube guarantees leakproof operation.
- The angled neck and the tip assure a uniform and extra fine stream of liquid.
- Unique dispensing tube design optimise backflow and turbulence free liquid transfer.

catalogue number	volume	cap colour	pack quantity
062.09.05R	500 ml	red	1 piece
062.09.05Y	500 ml	yellow	1 piece
062.09.05G	500 ml	green	1 piece
062.09.05B	500 ml	blue	1 piece

**"integral"**

- Manufactured from a mixture of high and low density polyethylene.
- Semi rigid, translucent and with high clarity with ergonomic shape for easy grip and handling.
- Single piece cap and delivery tube assembly eliminate air or liquid absorption through the cap and ensures a positive leaktight seal.
- Dispensing tip provides extra fine stream. The tip of the delivery tube can be cut in order to increase the liquid flow rate.
- Close the tip by means of the cap attached to the delivery tube.

catalogue number	volume	pack quantity
062.08.250	250 ml	1 piece
062.08.500	500 ml	1 piece



Dispensing tip provides extra fine stream. The tip of the delivery tube can be cut in order to increase the liquid flow rate.



Colour coded caps enables easy identification of the bottle content. Ribbed cap design allows easy opening and closing.



Dispensing tip provides extra fine stream. The tip of the delivery tube can be cut in order to increase the liquid flow rate.