



CARBOYS - "glass"

- Manufactured from borosilicate glass resistant to all chemicals according to ISO 4796.
- Ideal for safe storage, transportation and dispensing of liquids.
- Practical using thanks to wide GL 45 size mouth for easy filling.
- Supplied with a leak proof stopcock which prevents any vacuum formation when draining liquid out.

catalogue number	volume	quantity in pack
062.27.905	5 L	1 piece
062.27.910	10 L	1 piece
062.27.920	20 L	1 piece



CARBOYS - "P.P"

- Manufactured from polypropylene with UV stabilized heavy walls with stable carrying handles.
- Ideal for safe storage, transportation and dispensing of liquids.
- Wide neck for easy filling, and emptying. Featured with carrying handle.
- Supplied with a leak proof stopcock which prevents any vacuum formation when draining liquid out.

catalogue number	volume	quantity in pack
062.26.910	10 L	1 piece
062.26.920	20 L	1 piece



ASPIRATOR BOTTLES

- Manufactured from polyethylene with UV stabilized heavy walls with stable carrying handles.
- Ideal for safe storage, transportation and dispensing of liquids.
- Wide neck for easy filling, and emptying. Featured with carrying handle.
- Supplied with a leak proof stopcock which prevents any vacuum formation when draining liquid out.

catalogue number	volume	quantity in pack
062.28.905	5 L	1 piece
062.28.910	10 L	1 piece
062.28.920	20 L	1 piece



JERRYCANS - "polyethylene" - "with & W/O stopcock"

- Manufactured from blow molded robust HD-PE with heavy walls and handles. Provides the perfect solution for reliable and secured packaging, storing and distribution during transportation.
- Designed with good stackability features as well as leak-proof airtight which ensure non contamination of content to the users.



Cap with self sealing lock with integral ring.



Leak proof stopcock for draining liquid out.

catalogue number	volume	stopcock	base dimensions	height	pack quantity
062.21.902	2 L	no	150 x 82 mm	235 mm	1 piece
062.21.903	3 L	no	190 x 95 mm	270 mm	1 piece
062.21.905	5 L	no	183 x 138 mm	278 mm	1 piece
062.21.910	10 L	no	226 x 190 mm	325 mm	1 piece
062.21.920	20 L	yes	292 x 245 mm	359 mm	1 piece
062.21.930	30 L	yes	355 x 288 mm	395 mm	1 piece