



NEW

JARS - "wide neck" - "borosilicate glass"

- Manufactured from borosilicate glass with wide neck and volume inscriptions printed on side wall.
- Ideal for storing and transportation of both liquid and powder products. Supplied with polypropylene cap.

catalogue number	jar volume	base dimensions	jar height	pack quantity
061.03.500	500 ml	94 x 94 mm	141 mm	1 piece
061.03.901	1000 ml	105 x 105 mm	187 mm	1 piece
061.03.902	2000 ml	115 x 115 mm	260 mm	1 piece



NEW

JARS - "wide neck" - "neutral glass"

- Manufactured from neutral glass with wide neck and non corrosive metal cap. Ideal for storing and transportation of both liquid and powder products.

catalogue number	jar volume	neck width	jar height	pack quantity
061.29.040	40 ml	44 mm	50 mm	10 pieces
061.29.105	105 ml	49 mm	84 mm	10 pieces
061.29.370	370 ml	72 mm	121 mm	10 pieces
061.29.720	720 ml	90 mm	141 mm	10 pieces
061.29.901	1000 ml	94 mm	180 mm	10 pieces



NEW

JARS - "wide neck" - "P.P"

- Manufactured from clear polypropylene with wide neck and polypropylene screw cap with silicone ring.
- Ideal for storing of both liquid and powder products.

catalogue number	jar volume	dimensions Ø x height	pack quantity
061.13.001	250 ml	90 x 62 mm	1 piece
061.13.002	410 ml	90 x 102 mm	1 piece
061.13.004	1270 ml	140 x 110 mm	1 piece
061.13.006	1420 ml	140 x 135 mm	1 piece
061.13.008	1800 ml	140 x 160 mm	1 piece



NEW

BOTTLES - "wide neck" - "neutral glass"

- Manufactured from neutral glass with round neck according to DIN 168. Non autoclavable. Supplied with cap.
- Ideal for storing and transportation of both liquid and powders. Supplied with polypropylene cap with insert.
- Amber bottles prevent from UV lights to penetrate inside the bottle.

bottles clear	bottles amber	bottle volume	bottle diameter	bottle height	quantity in pack
061.10.050	061.11.050	50 ml	44 mm	75 mm	10 pieces
061.10.100	061.11.100	100 ml	50 mm	93 mm	10 pieces
061.10.250	061.11.250	250 ml	70 mm	113 mm	10 pieces
061.10.500	061.11.500	500 ml	85 mm	154 mm	10 pieces
061.10.901	061.11.901	1000 ml	104 mm	181 mm	10 pieces