



EVAPORATING DISHES - "porcelain" - "flat"

➤ Manufactured from porcelain as flat bottom according to DIN 12907. Externally and internally manufactured from glazed porcelain except the spout.

➤ A perfect product for heating and evaporating liquids.

➤ Not recommended to use with direct flame in order to prevent breakage.

catalogue number	external diameter	height	volume	pack quantity
038.11.060	60 mm	13 mm	20 ml	6 pieces
038.11.080	80 mm	16 mm	60 ml	6 pieces
038.11.100	100 mm	20 mm	125 ml	6 pieces



EVAPORATING DISHES - "porcelain" - "round"

➤ Manufactured from porcelain as round bottom according to DIN 12907. Externally and internally manufactured from glazed porcelain except the spout.

➤ A perfect product for heating and evaporating liquids. Not for use with direct flame in order to prevent breakage.

catalogue number	external diameter	height	volume	pack quantity
038.12.060	60 mm	25 mm	30 ml	6 pieces
038.12.080	80 mm	32 mm	80 ml	6 pieces
038.12.100	100 mm	40 mm	175 ml	6 pieces
038.12.125	125 mm	50 mm	450 ml	6 pieces
038.12.150	150 mm	65 mm	750 ml	6 pieces



EVAPORATING DISHES - "st. steel" - "round"

➤ Manufactured from stainless steel as round bottom.

➤ A perfect product for heating and evaporating liquids.

➤ Can be used with direct flame in order to prevent breakage.

catalogue number	external diameter	height	volume	pack quantity
038.14.060	60 mm	30 mm	50 ml	1 piece
038.14.080	80 mm	40 mm	130 ml	1 piece
038.14.100	100 mm	50 mm	250 ml	1 piece



EVAPORATING DISHES - "quartz glass"

➤ Manufactured from quartz glass as round bottom.

➤ A perfect product for heating and evaporating liquids.

➤ Can be used with direct flame in order to prevent breakage.

catalogue number	external diameter	height	volume	pack quantity
038.16.060	60 mm	32 mm	20 ml	1 piece
038.16.080	80 mm	45 mm	49 ml	1 piece
038.16.100	100 mm	55 mm	90 ml	1 piece