



### CRYSTALLIZATION DISHES - "glass"

➤ Manufactured from borosilicate glass 3.3 according to DIN 12338.

➤ Used for growing crystals and any other application requiring a dish which is in height between petri dishes and glass beakers.

catalogue number	external diameter	height	pack quantity
049.05.060	60 mm	35 mm	1 piece
049.05.100	70 mm	40 mm	1 piece
049.05.150	80 mm	45 mm	1 piece
049.05.300	95 mm	55 mm	1 piece
049.05.900	140 mm	75 mm	1 piece



### CASSEROLES

➤ Manufactured from glazed porcelain according to DIN 12907.

catalogue number	external diameter	height	pack quantity
038.10.080	80 ± 10 mm	45 mm	1 piece
038.10.100	100 ± 10 mm	65 mm	1 piece



### CRYSTALLIZATION DISHES - "stainless steel"

➤ Alternative to borosilicate glass dishes for similar applications.

➤ Manufactured from seamless and robust stainless steel.

catalogue number	external diameter	height	pack quantity
049.08.070	70 mm	35 mm	1 piece
049.08.085	85 mm	45 mm	1 piece
049.08.110	110 mm	55 mm	1 piece



### COMBUSTION BOATS

➤ Manufactured from porcelain with temperature resistivity till 1050°C.

catalogue number	dimensions	pack quantity
038.17.001	80 x 13 x 8 mm	10 pieces
038.17.002	102 x 20 x 13 mm	10 pieces



### DISH TONG

➤ Manufactured from polished stainless steel.

➤ Due to unique design can hold all size of evaporating dishes.

➤ Curved teeth of the dish tong enables easy and safe handling of the evaporating dishes.

catalogue number	total length	pack quantity
038.05.001	280 ± 10 mm	1 piece