



CRUCIBLES - "porcelain" - "1050°C"

- Manufactured from glazed porcelain both inside and outside according to DIN 12904 and ISO 1772 standards.
- Porcelain crucibles can resist temperatures up to 1050°C.
- Supplied without cover lids.
- Cover lids should be ordered separately.

catalogue number	external diameter	height	volume	pack quantity
038.13.H30	30 mm	38 mm	15 ml	6 pieces
038.13.M35	35 mm	28 mm	13 ml	6 pieces
038.13.H35	35 mm	44 mm	24 ml	6 pieces
038.13.L40	40 mm	25 mm	15 ml	6 pieces
038.13.M40	40 mm	32 mm	20 ml	6 pieces
038.13.H40	40 mm	50 mm	35 ml	6 pieces
038.13.L45	45 mm	28 mm	25 ml	6 pieces
038.13.M45	45 mm	36 mm	30 ml	6 pieces
038.13.H45	45 mm	56 mm	50 ml	6 pieces
038.13.L50	50 mm	32 mm	29 ml	6 pieces
038.13.M50	50 mm	40 mm	39 ml	6 pieces
038.13.H50	50 mm	65 mm	75 ml	6 pieces
038.13.L60	60 mm	38 mm	45 ml	6 pieces
038.13.M60	60 mm	50 mm	85 ml	6 pieces
038.13.M70	70 mm	60 mm	140 ml	6 pieces

CRUCIBLES - "porcelain" - "1150°C"

- HALDENWANGER crucibles manufactured from glazed porcelain both inside and outside according to DIN 12904 and ISO 1772 standards.
- Due to the inherent physical characteristics crucibles can resist temperatures up to 1150°C. Supplied without cover lids.

catalogue number	manufacturer cat. no.	external diameter	height	volume	pack quantity
038.15.790	79C/0	30 mm	38 mm	15 ml	5 pieces
038.15.791	79C/1	35 mm	44 mm	24 ml	5 pieces
038.15.792	79C/2	40 mm	50 mm	40 ml	5 pieces



CRUCIBLE LIDS

- Manufactured from glazed porcelain.
- Porcelain crucible lids can resist temperatures up to 1050°C.

catalogue number	internal diameter	pack quantity
038.15.030	30 mm	6 pieces
038.15.035	35 mm	6 pieces
038.15.040	40 mm	6 pieces
038.15.045	45 mm	6 pieces
038.15.050	50 mm	6 pieces
038.15.060	60 mm	6 pieces

