



NEW

CRUCIBLES - "alumina"

➤ Alumina crucibles can resist to 1700°C.

catalogue number	external dia. x height	pack quantity
038.13.A20	33 x 44 mm	6 pieces
038.13.A50	46 x 53 mm	6 pieces
038.13.A99	64 x 86 mm	6 pieces



NEW

CRUCIBLES - "quartz"

➤ Quartz crucibles can resist to 1050°C.

catalogue number	external dia. x height	pack quantity
038.13.Q15	38 x 30 mm	1 piece
038.13.Q24	44 x 35 mm	1 piece
038.13.Q40	50 x 40 mm	1 piece



NEW

CRUCIBLES - "st. steel"

➤ Stainless steel crucibles can resist to 1200°C.

catalogue number	external dia. x height	pack quantity
038.13.S35	40 x 35 mm	1 piece
038.13.S45	50 x 45 mm	1 piece
038.13.S50	60 x 50 mm	1 piece

TONG - "crucible"

- Manufactured from stainless steel.
- Curved teeth of the tong provides easy handling.

catalogue number	total length	pack quantity
038.06.250	250 ± 10 mm	1 piece
038.06.300	300 ± 10 mm	1 piece
038.06.400	400 ± 10 mm	1 piece

