

BUNSEN BURNERS - "for all gases such as neutral gas, bottle gas, coal gas"



"mini"

- Low form standard Bunsen burner for standard flame.
- Gas flow controlled by a rotatable regulator for easy flame adjustment.
- With built-in flame retainer and adjustable air regulator.
- Mini design allows usage even in short and narrow places.

catalogue number	pack quantity
050.01.005	1 piece



"standard"

- Standard Bunsen burner for standard flame.
- Gas flow controlled by a rotatable regulator for easy flame adjustment.
- With built-in flame retainer and adjustable air regulator.
- The stabilizer top allows gas to emerge at a low velocity in an annular ring acting as a pilot for the regular flame.

catalogue number	pack quantity
050.01.001	1 piece



"tirril"

- Tirril burners allow greater flexibility in the adjustment of the air-gas mixture.
- Easy adjustment in the air&gas mixture even for very slow heating or quick boiling.
- Suitable for high temperatures up to 1500°C and with its compact design, can be used with most of the tripods.

catalogue number	pack quantity
050.01.002	1 piece



"meker"

- Suitable for high temperature work up to 1700°C.
- Has a perforated plate of 25 mm diameter which provides flame stability and ensures high and evenly distributed flame.
- Gas flow controlled by a rotatable regulator for easy flame adjustment.
- Other technical specifications are same as other Bunsen burners.

catalogue number	pack quantity
050.01.003	1 piece



BUNSEN BURNER HOSE

- Manufactured from thick rubber hose specially designed to be used with Bunsen burners. Reinforced ends fit to Bunsen burner and gas tap securely without any risk of gas leakage.

catalogue number	hose wall thickness	hose length	pack quantity
050.08.075	6 mm	75 cm	1 piece
050.08.100	6 mm	100 cm	1 piece



BUNSEN BURNER STAND

- Manufactured from spray coated heavy casted metal.
- Perfectly fits to burner base and prevents the burner tilt down due to the torsion and twisting effect of burner hose.
- Perfectly fits to the base shape of Bunsen burners 050.01.002 and 050.01.003.

catalogue number	pack quantity
050.09.001	1 piece