

CLAMPS - "2 prongs"



- Manufactured from chrome plated cast steel with two prongs which are PVC coated for holding the objects gently.
- Two separately adjustable clamps enables perfect control of holding pressure and precise positioning of the objects.
- Clamps can be opened and closed smoothly with a grip width of 0-50 mm.

<i>solid stem clamp</i>	<i>swivel clamp</i>	<i>jaw opening</i>	<i>pack quantity</i>
045.01.001	045.04.001	0 - 50 mm	1 piece

- Available 2 types:
 - with solid stem of 12 mm rod diameter with the holding angle and distance from the support rod can be adjusted by a separate bosshead.
 - with swivel mechanism which enables the object freely swivel 360° and lockable at any angle. Complete with a integral bosshead usable for rods maximum diameter of 21 mm.

CLAMPS - "3 prongs"



- Manufactured from chrome plated cast steel with two prongs which are PVC coated for gentle hold of the objects.
- Two separately adjustable clamps enables perfect control of holding pressure and precise positioning of the objects.
- Clamps can be opened and closed smoothly with a grip width of 0-70 mm.

<i>solid stem clamp</i>	<i>swivel clamp</i>	<i>jaw opening</i>	<i>pack quantity</i>
045.03.001	045.06.001	0 - 70 mm	1 piece

- Available 2 types:
 - with solid stem of 12 mm rod diameter with the holding angle and distance from the support rod can be adjusted by a separate bosshead.
 - with swivel mechanism which enables the object freely swivel 360° and lockable at any desired angle. Complete with an integral bosshead usable for rods max.