


**BOSSHEADS - "standard use"**

- Manufactured from chrome plated cast alloy with smoothly opening distance of 0 - 15 mm.
- With fix stem that holds the objects on fixed position.

catalogue number	suitable for rods	pack quantity
046.01.002	0 - 16 mm	1 piece

**BOSSHEADS - "standard use" - "swivel"**

- Manufactured from chrome plated cast alloy with smoothly opening distance of 0 - 15 mm.
- With swivel mechanism which enables the object freely swivel 360° and lockable at any desired angle.

catalogue number	suitable for rods	pack quantity
046.01.001	0 - 16 mm	1 piece

**BOSSHEADS - "heavy duty"**

- Small size bosshead manufactured from chrome plated cast steel with smoothly opening distance of 0 - 16 mm.
- Due to the sturdy construction specially suitable to hold heavy objects.

catalogue number	size	suitable for rods	pack quantity
046.05.001	small	0 - 16 mm	1 piece

**BOSSHEADS - "heavy duty"**

- Large size bosshead manufactured from powder coated cast steel with smoothly opening distance of 0 - 20 mm.
- Due to the sturdy construction specially suitable to hold heavy objects.

catalogue number	size	suitable for rods	pack quantity
046.05.002	large	0 - 20 mm	1 piece

**BOSSHEADS - "for frame work"**

- Manufactured from chrome plated cast steel and suitable for frame works and fixing perpendiculars.
- Can be moved vertically and horizontally along the rods.

catalogue number	suitable for rods	pack quantity
046.07.001	12 mm	1 piece

