

## SPATULAS - "macro spoon"



➤ Manufactured from high quality (AISI 304) polished stainless steel. Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

➤ Most commonly used spatula for general laboratory applications.

➤ One end flat spatula and the other end as macro spoon.

catalogue number	spatula total length	pack quantity
047.01.120	120 ± 5 mm	1 piece
047.01.150	150 ± 5 mm	1 piece
047.01.180	180 ± 5 mm	1 piece
047.01.210	210 ± 5 mm	1 piece
047.01.250	250 ± 10 mm	1 piece
047.01.300	300 ± 10 mm	1 piece

## SPATULAS - "chattaway"



➤ Manufactured from high quality (AISI 304) polished stainless steel. Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

➤ Specially efficient and practical to hold products from flat and sticky surfaces due to the bent flat end.

➤ One end flat and the other end as bent flat.

catalogue number	spatula total length	pack quantity
047.09.100	100 ± 5 mm	1 piece
047.09.130	130 ± 5 mm	1 piece
047.09.150	150 ± 5 mm	1 piece
047.09.180	180 ± 5 mm	1 piece
047.09.210	210 ± 5 mm	1 piece

## SPATULAS - "poly"



➤ Manufactured from high quality (AISI 304) polished stainless steel. Robust seamless construction guarantees long working life and smooth finish enables easy cleaning.

➤ Specially efficient and practical to hold products from the edges of flasks due to the ellipsoidal grooved end.

➤ One end flat spatula and the other end as grooved.

catalogue number	spatula total length	pack quantity
047.08.150	150 ± 5 mm	1 piece
047.08.180	180 ± 5 mm	1 piece
047.08.210	210 ± 5 mm	1 piece