



SCOOPS - "P.P" - "set of 3"

- Manufactured from polypropylene. Volume marked.
- Scoops have flat bottom provides to stand horizontally.
- Scoops have rounded corners for easy cleaning.
- Excellent for scooping of small or large volumes anything from powder to ice.
- Packed as 3 scoops as a set with scoop volume of 50, 100 and 250 ml.

catalogue number	pack quantity
037.26.001	1 set



SCOOPS - "P.E" - "macro"

- Manufactured from HD-PE.
- Scoops have rounded corners for easy cleaning.
- Excellent for scooping of small or large volumes anything from powder to ice.

catalogue number	size	scoop volume	pack quantity
037.27.001	small	500 ml	1 piece
037.27.002	medium	1000 ml	1 piece
037.27.003	large	2000 ml	1 piece

The indicated volumes are an estimated volume and are not meant for taking exact sample volumes.



SCOOPS - "stainless steel"

- Manufactured from stainless steel.
- Scoops have rounded corners for easy cleaning.
- Excellent for scooping of small or large volumes anything from powder to ice.

catalogue number	scoop diameter	pack quantity
037.08.050	50 ml	1 piece
037.08.075	75 ml	1 piece
037.08.100	100 ml	1 piece
037.08.120	120 ml	1 piece

The indicated volumes are an estimated volume and are not meant for taking exact sample volumes.