

single channel micropipettes

non autoclavable



- ▶ High quality micro pipettes with reasonable prices to cover complete pipetting range.
- ▶ It is constructed of high quality plastic and combined with stainless steel piston, thus the item is robust yet lightweight.
- ▶ Individually calibrated according to **ISO 8655 standards** and the quality certificate for accuracy and precision is provided with each pipette.
- ▶ The ergonomical, rounded hand grip and finger handle ensure comfortability and stability as well as a correct pipetting angle.
- ▶ The unique design allows quick and easy volume setting with only one hand.
- ▶ Built-in streamlined tip ejector enables easy tip ejection and allows access to narrow neck volumetric flasks and to the bottom of the tubes.
- ▶ Click stop digital system for volume setting. The volume display is clearly visible all the time during operation.

non autoclavable pipettes	volume range	test volume	Inacc. %	Impr. %
001.01.002	0.2 - 2.0 µl	0.2 µl	12.0 %	6.0 %
		1 µl	2.5 %	1.5 %
		2 µl	1.4 %	1.7 %
001.01.010	0.5 - 10 µl	0.5 µl	5.0 %	2.8 %
		5 µl	2.5 %	0.8 %
		10 µl	1.5 %	0.4 %
001.01.020	2 - 20 µl	2 µl	6.0 %	5.0 %
		10 µl	1.0 %	0.6 %
		20 µl	0.8 %	0.3 %
001.01.050	5 - 50 µl	5 µl	2.0 %	0.5 %
		25 µl	1.2 %	0.3 %
		50 µl	0.7 %	0.3 %
001.01.100	10 - 100 µl	10 µl	2.0 %	0.3 %
		50 µl	0.7 %	0.3 %
		100 µl	0.8 %	0.2 %
001.01.200	20 - 200 µl	20 µl	1.5 %	0.3 %
		100 µl	0.4 %	0.2 %
		200 µl	0.8 %	0.2 %
001.01.901	100 - 1000 µl	100 µl	0.6 %	0.3 %
		500 µl	0.3 %	0.2 %
		1000 µl	0.5 %	0.2 %
001.01.905	1000 - 5000 µl	1000 µl	0.8 %	0.2 %
		2500 µl	0.3 %	0.2 %
		5000 µl	0.6 %	0.2 %
001.01.910	1000 - 10000 µl	1000 µl	0.8 %	0.2 %
		5000 µl	0.3 %	0.2 %
		10000 µl	0.6 %	0.2 %