

wide neck bottles polyethylene




- Manufactured from high density polyethylene which can resist high pressures before rupture.
- These space saving bottles are ideal for long term storage of products.
- Approximate graduations allow filling without the usage of cylinders.
- Amber bottles are manufactured from a special resin which prevents UV lights penetrate inside the bottle.
- Inscriptions and graduations are in embossed form.



The head space between cap and insert plug can be used for private and confidential notes.

SECURITY SEAL



Security seal is used ensure the safety locking and preserveness of the product inside the bottles. 4 different attractive colors are used to identify among the seals which have unique identification numbers printed on the seal plate.



clear bottles	amber bottles	volume	clear in box	amber in box
061.07.050	061.08.050	50 ml	10 pcs	1 piece
061.07.100	061.08.100	100 ml	10 pcs	1 piece
061.07.250	061.08.250	250 ml	10 pcs	1 piece
061.07.500	061.08.500	500 ml	10 pcs	1 piece
061.07.901	061.08.901	1000 ml	10 pcs	1 piece
061.07.902	061.08.902	2000 ml	6 pcs	1 piece



Insert the stainless steel seal wire through the holes of both the cap and bottle itself. Then insert the wire into the hole of the seal as shown above.

Last step is to wind the seal key clockwise until the seal wire gets tight enough to provide optimum security lock.
Catalogue No: 061.09.001