

thermometers

enclosed scale

- Enclosed Scale Thermometers are double walled thermometers with mercury column.
- Enclosed white scale enables accurate reading of temperature.

catalogue number	operation range °C	subdivision	quantity in pack
059.01.001	-10 / + 50	1°C	1 piece
059.01.002	-10 / +110	1°C	1 piece
059.01.003	-10 / +150	1°C	1 piece
059.01.004	-10 / +200	1°C	1 piece
059.01.005	-10 / +250	1°C	1 piece
059.01.006	-10 / +300	2°C	1 piece
059.01.007	-10 / +360	2°C	1 piece
059.01.008	-10 / +420	2°C	1 piece

thermometers

solid stem

- Solid Stem Glass Thermometers are calibrated for partial immersion depth of 76 mm.
- Thermometer has red alcohol column and the body is reflecting white colour which provides easier readings.

catalogue number	operation range °C	subdivision	quantity in pack
0059.02.001	-20 / + 50	0.5 °C	10 pieces
059.02.002	-20 / +110	0.5 °C	10 pieces
059.02.003	-10 / +150	1 °C	10 pieces
059.02.004	-10 / +200	1 °C	10 pieces

thermometer clamp

with bosshead

- Thermometer Clamp & Holder with boss head has an opening of 0-30 mm.
- Steel spring arm holds the thermometer firmly and keeps the thermometer at the center.
- With swivel mechanism which enables the object to swivel 360° and is lockable at any desired angle.

Catalogue No: 060.04.001