

NSZ-806 Zoom Stereo Microscope



Description

- Excellent Image Quality & Optical Performance with Infinity Parallel Optical System.
- Comfortable Operation with Ergonomic Design Principle.
- With LED light for both incident and transmitted illumination, Providing even illumination and life expectancy can reach 6000 hours.
- A Wide Range of Accessories to be Equipped with the Microscope for Multi-purposelt is an Outstanding Assistant in Industrial, Medical and Scientific Research Areas.

Microteb

Tel: 021 66 97 60 37

Fax: 021 66 90 33 54

www.microteb.com

Specification:

		NSZ-806
Drawing Tube	Binocular Head, 20°Inclination	●
	Tilting Binocular Head, 0°-35° Inclination	
	Photo Attachment(one photo tube)	
	Photo Attachment(two photo tube)	○
Eyepiece	EW10×/22	●
	WF15×/16	○
	WF20×/12	○
Zoom Objective	0.8×-5×	●
	0.8×-6.4×	
	0.8×-8×	
Objective	Plan Achromatic Objective 1×	●
	Plan Apochromatic Objective 0.5×	○
	Plan Achromatic Objective 2×	○
	Achromatic Objective 0.3×	
Zoom Ratio	1: 6	●
	1: 8	
	1: 10	
Working Distance	78mm	●
Focusing	Focusing Range 105mm	●
Illumination	Transmission or Reflection LED Illumination, Brightness Adjustable	●
	Coaxial Illumination	○
Fluorescent Attachment		○
Dark Field Attachment		○
Photo/ Video Attachment		○
Simple Polarization Set		○
Iris Diaphragm		○
Mechanical Stage		○
Temperature Stage		○
Universal Table Stand		○
D Stand		○
Coaxial Coarse & Fine Focusing Unit		○

Note: ● Standard Outfit, ○ Optional